

CLASSIFIED ADVERTISEMENTS

The rate for classified advertisements is 5s. a line (or space equivalent of a line), with an extra charge of 1s. for the use of a Box Number. Semi-displayed classified advertisements are 60s. per single-column inch.

CHEMIST—required in Analytical Department of an organisation with wide interests, particularly in food stuffs. The work is varied and non-routine, and experience in a Public Analyst's or similar laboratory would be considered an advantage. Applicants should be of graduate standard, but sound analytical technique is more important than academic qualifications. Five-day week, non-contributory Pension Scheme. Write, giving age and particulars of qualifications and experience to Box 4048, The Analyst, 47 Gresham Street, London, E.C.2.

ANALYTICAL CHEMIST

(Honours Degree, Grad. R.I.C. or A.R.I.C.)

A vacancy occurs in the Analytical Laboratory of the Research and Development Department of BAKELITE LIMITED for a qualified Analytical Chemist. The work is interesting and non-routine and includes analysis of a wide variety of plastics, their raw materials and associated products.

Excellent working conditions in a progressive organisation as well as superannuation benefits and social and welfare activities.

Apply to Research and Development Manager, Bakelite Limited, Tyseley, Birmingham, 11, giving age, qualifications and experience.

VACANCY in Consulting Laboratories in N.W. London for Senior Chemist with experience in General Food Analysis. Apply Box No. 4050, The Analyst, 47 Gresham Street, London, E.C.2.

FREDK. BOEHM LIMITED, require a **Development Chemist**, to work on problems involving Protein Hydrolysates and Synthetic Resins and an **Analyst** to take charge of Quality Control Laboratory. Applicants should be aged 22-30. Experience in this particular type of work is desirable but at least two years' industrial experience is essential.

Qualifications, A.R.I.C., grad. R.I.C. or B.Sc.(Hons). In the case of the analytical chemist a knowledge of Bacteriology and Microbiological assay technique would be advantageous. Good salary and prospects.

Applications should be addressed to:—

Chief Food Chemist,
Fredk. Boehm Ltd.,
Crabtree Manorway,
Belvedere, Kent.

SCHWEPPE'S (HOME) LIMITED, require a Junior Chemist to work in their Production Laboratory. Applicants should have O.N.C. or at least "A" level in Chemistry. Hours of work 7.30 a.m. to 4.45 p.m., Monday to Thursday, with a 4.15 p.m. finish on Friday. During certain periods applicants might be expected to work limited shift. (Not nights) 5-day week. Sick Pay, Pension and Profit Sharing schemes in operation.

Please apply in writing to Personnel Manager, Schweppe's (Home) Ltd., Cray Road, Sidcup.

ANALYSTS REQUIRED with experience of foods and drugs, preferably in Public Analysts' laboratories. Applicants should have at least HNC or equivalent qualifications. A wide range of analytical work demands competence and versatility but offers good prospects in a rapidly expanding organisation. Opportunities to study for higher qualifications. Apply: J. A. Radley (Laboratories) Ltd., 220 Elgar Road, Reading, Berks.

KENT COUNTY COUNCIL require **ASSISTANT ANALYST** in the County Laboratory. Graduate or equivalent. Salary scale £960-£1,140. Further particulars from County Analyst, County Hall, Maidstone.

CITY AND COUNTY OF KINGSTON UPON HULL**DEPUTY CITY ANALYST**

Applications are invited for the appointment of Deputy City Analyst, the salary to be within J.N.C. Grade A, viz. £1,365 rising to £1,565 per annum.

Particulars of the appointment and forms on which application should be made are obtainable from the City Analyst, 184 High Street, Kingston upon Hull, by whom applications must be received not later than April 3rd, 1962.

IMPERIAL COLLEGE OF SCIENCE AND TECHNOLOGY

A Bursary, value not less than £400 p.a. plus fees, suitable for graduate in Chemistry or Metallurgy who wishes to take higher degree, will become available in October, 1962. The research will be on the polarography of inorganic cations in organic solvents with a view to application in the analysis of ferrous metallurgical materials. Applications, with summary of career and names of referees, should be sent to Dr. D. A. Pantony, Metallurgy Dept., Imperial College, London, S.W.7.

RIKER LABORATORIES**Pharmaceutical Manufacturers**

Have a vacancy for an

ANALYTICAL CHEMIST

The scope of the work is wide and interesting and includes all analytical aspects relating to raw materials, intermediates and finished pharmaceutical products.

Applicants should be under 35 years of age and possess a B.Sc., B.Pharm. or A.R.I.C. qualification. Some previous experience in the pharmaceutical industry would be an advantage and familiarity with modern instrumental techniques is desirable. Applications from recently qualified analysts will also be considered.

The Company is progressive and rapidly expanding and provides excellent working facilities. There is a contributory pension scheme, and salary will be commensurate with age, qualifications and experience.

Apply in writing to the Personnel Director, Riker Laboratories, Morley Street, Loughborough, Leics.

A QUALIFIED CHEMIST is required to fill a vacancy in the Chemical Laboratory. Experience in general analytical techniques and a bias towards organic chemistry is desirable. Starting salary £700-£800 p.a. Applications naming two referees should be addressed to the Director, Fruit and Vegetable Canning and Quick Freezing Research Association, Chipping Campden, Glos.

EXPERIENCED ANALYST

Johnson, Matthey & Co., Limited, chemists and metallurgists specialising in the precious and rare metals, have a vacancy for a qualified experienced analytical chemist of degree standard.

The position is in the Company's laboratories at a new chemical production centre at Royston, North Hertfordshire and some 12 miles from Cambridge. The laboratories include Spectrographic and X-Ray Fluorescence departments.

The successful candidate will be required to take charge of a chemical laboratory, to carry out referee analysis, and to develop new methods of analysis, in a very wide and interesting range of precious metals, ores, catalysts, rare earths and many other rare elements. The Company is not, however, expecting to find specialised experience in these particular fields.

Applications should be made in writing, giving full details of qualifications, experience, age and an indication of present position and salary, to the Company Secretary, 78 Hatton Garden, London, E.C.1, quoting reference B.31.



ANALYTICAL CHEMIST

Applications are invited, from men or women, for a position as Analytical Chemist in the laboratories at Bournville.

Applicants should be qualified, B.Sc., G.R.I.C., or A.R.I.C. Preference would be given to those with experience in the food industry or a Public Analyst's laboratory.

The successful applicant will be involved almost entirely in analytical development work and the investigation of new techniques.

The salary offered will match qualifications and experience. Conditions of employment are excellent and there is a good Pension Scheme.

Applications should be addressed to:—

**THE CHIEF CHEMIST,
CADBURY BROS. LTD.,
BOURNVILLE, BIRMINGHAM**

DELANEY GALLEY require an Analytical Chemist preferably A.R.I.C., B.Sc. or Dip. Tech. for their Cricklewood Laboratory. The Company specialises in heat exchange and metal joining. Ideally the successful applicant will be under 29 years and a British National. Please apply in writing to Personnel Officer, Delaney Galley Ltd., Vulcan Works, Edgware Road, Cricklewood, N.W.2.

EDITORIAL ASSISTANT required in the office of *Analytical Abstracts*. Some knowledge of chemistry is desirable. Five-day week. Hours 9.30 to 5. Apply in writing to The Editor, *Analytical Abstracts*, 14 Belgrave Square, London, S.W.1.

**BRITISH BAKING INDUSTRIES RESEARCH
ASSOCIATION**
require

Scientific Officer, male, graduate or A.R.I.C. for Analytical Laboratory. Duties will include research into analytical methods and the investigation of a wide range of problems associated with the baking industry. Experience of food analysis an advantage. This is a responsible post in a vigorous organisation. Pensionable under F.S.S.U.

Applications should include details of age, qualifications, experience and names of two referees. Write to:

The Director of Research, British Baking Industries Research Association, Chorleywood, Rickmansworth, Herts.

PUBLIC ANALYST requiring **DEPUTY** invites applications for post of Senior Assistant from persons having the necessary qualifications. Branch E certificate essential. Salary will be commensurate with qualifications and experience and the post offers excellent prospects for quick advancement.

Applications to Dr. H. K. Lawton, F.R.I.C., Room 74, 16 Donegall Square South, Belfast, Northern Ireland, from whom further particulars may be obtained.

**BRITISH CHEMICAL
&
SPECTROGRAPHIC
STANDARDS**

see our
EXHIBITS
STAND 51

**ANALOID REAGENTS
FOR
METALLURGICAL
ANALYSIS**

LABORATORY APPARATUS & MATERIALS EXHIBITION

14th — 18th MAY

THE EXHIBITION HALL, HARROGATE, YORKS.

**BUREAU OF ANALYSED SAMPLES LTD. I RIDSDALE & CO. LTD.
NEWHAM HALL, MIDDLESBROUGH**



MULLARD
RESEARCH LABORATORIES
Solid State Physics Division
ANALYST/SPECTROGRAPHER

required to establish a new laboratory for the spectrographic analysis of a variety of new and unusual materials used in research on semiconductors, magnetic oxides, master materials, etc.


The post presents an opportunity for a man under thirty with an honours degree or its equivalent who is able to apply a sound knowledge of physics and chemistry to the development of new techniques. Some experience of emission spectrography or other physical methods of analysis is desirable, and the work will require a critical approach to the assessment of available methods for the determination of impurities at low level and an inventive approach to instrumentation.

Applications, ref. KH/1 to:

Mr. G. A. Taylor
 Mullard Research Laboratories,
 Salfords,
 nr. Redhill, Surrey.

perHaps?

**NO! YOU CAN BE CERTAIN OF
 ACCURATE PH READINGS WITH**



JOHNSON
INDICATOR PAPERS

Johnson COMPARATOR TEST PAPERS give a wide range of readings from 3.6 to 10 pH in steps of 0.3 pH and from 1.0 to 3.5 in steps of 0.5 pH. JOHNSON UNIVERSAL TEST PAPERS give the pH value of any solution to within 0.5 pH in the range of 1 to 11. These papers have long been accepted by Science and Industry for their consistent accuracy. Johnsons also supply LITMUS paper and a series of nearly forty other Indicator Papers for many specific test purposes.

Write for leaflet

JOHNSONS OF HENDON LTD.

LONDON N.W.4

BINDING

Have your back numbers of
The Analyst bound in the
 standard binding case.
 Send the parts together with
 a remittance for 16/6

To:—

W. HEFFER & SONS LTD
HILLS ROAD CAMBRIDGE

DECENNIAL INDEXES TO THE ANALYST

1926-1935	(Sewn) £1 5s. net
1936-1945	(Cloth) £2 2s. net
1946-1955	(Cloth) £5 net

*It is regretted that earlier volumes are
now out of print.*

(Available at reduced prices to members of the Society for
Analytical Chemistry on direct application to The Editor,
The Analyst, 14 Belgrave Square, London, S.W.1)

Published for the Society by

W. HEFFER & SONS LTD., CAMBRIDGE

Reprint of **THE ANALYST**


VOLUME I (1876) to VOLUME X (1885)

Bound in cloth

Each Volume **£6 6s.** net

The early volumes of this indispensable journal, out of print soon
after publication and lacking from many library sets, have now
been made available again in response to repeated requests.
We hope to reprint all unavailable volumes, if the demand
proves adequate.

W. HEFFER & SONS, LTD.
CAMBRIDGE — ENGLAND



NEW WAY TO HUNT THE METALLURGICAL THIMBLE



PHILIPS NEW GAMMA SPECTROMETER

Which metallurgical thimbles do you hunt?

Material distribution . . . structural formations . . . deformations? For these, and similar aspects of metallurgical application analysis, you now have a new and valuable additional tool: Philips Gamma Spectrometer. Highly sensitive, consistently reliable, Philips new Gamma Spectrometer is designed to give accurate evaluations of materials labelled with radio-active isotopes. It will trace radio-active source materials in sheet, bar or billet metals and alloys without destroying them.

WIDER RANGE—AUTOMATIC OPERATION. The new Philips Gamma Spectrometer incorporates a single channel analyser providing bias and discriminator levels up to 100 volts, channel widths up to 32 volts.

It also uses the Philips Non-overloading Linear Amplifier. Automatic scanning of the bias range, in forward or reverse directions at four alternative speeds, and automatic indications on completion of scan, coupled with recorder control facilities, allow the operator to leave the equipment unattended. The high-stability EHT unit provides negative or positive supplies. Gain Scanning is also provided for in this equipment.

This new Philips Spectrometer gives the standards required by research, with the reliability and reasonable price demanded by industry. It is built up from the Philips Unitised range of Nucleonic instruments, and users of these units can assemble the Spectrometer by appropriate addition. Write for full details.

FREE APPLICATION SERVICE. To check the value of this equipment in relation to your own problems, send a sample for preliminary analysis by our specialists.

Overseas enquiries, please, direct to the manufacturers: N.V. Philips, Eindhoven, The Netherlands.

Sole distributors in U.K.:

RESEARCH & CONTROL INSTRUMENTS LTD.

Instrument House, 207 King's Cross Road, London, W.C.1. Telephone TERminus 2877 (RCL0060A)

*more & more tubing
in production!*

VITREOSIL
Registered Trade Mark

pure fused quartz

New production techniques are now speeding up deliveries of Vitreosil Fused Quartz Tubing—the most transparent manufactured material. It is impermeable, can be used repeatedly up to 1050°C, is inert to acids except hydrofluoric, and has high light transmittance in the U.V. and I.R. spectra.

all BRITISH



...the BEST

There are a wide range of diameters and two wall thicknesses—thin and normal.

Full details on request.

THE THERMAL SYNDICATE LTD.

P.O. Box No. 6, Wallsend, Northumberland.

Tel: Wallsend 62-3242/3

LONDON: 9, BERKELEY STREET, W.1.


Tel: Hyde Park 1711/2



A Major Advance in Moisture Determination

M&B Karl Fischer Reagent **stabilized**

Titre 4 mg. of water per ml.



- Single solution in 500 ml. bottles
- About fifteen times more stable than the conventional reagent
- Exceptional stability extends scope to the occasional user
- Restandardization once a fortnight instead of once a day
- More than 90 per cent time and reagent saved on restandardizations
- Solution in the reservoir remains serviceable for months


Full details are available on request

M&B
brand laboratory chemicals and reagents

Manufactured by **MAY & BAKER LTD**

DAGENHAM · ESSEX · TEL: DOMINION 3060 EXTS: 320 and 316

LA694



Matthey

GRAPHITE

FOR SPECTROGRAPHY

Now available in the Matthey range of Spectrographically Standardised Substances, is a new grade of highly-purified graphite.

Of excellent physical and chemical properties, the new graphite is available in the form of rod and conforms to the following specification for purity:

Total impurities detected spectrographically do not exceed 5 parts per million.

Analyses of current production show total impurities of about 2 parts per million only.

A second grade of graphite, specially prepared for use in analyses of a routine nature, conforms to the following:


Total impurities detected spectrographically do not exceed 100 parts per million.


On typical batches total impurities are in the range 30-60 parts per million and consist primarily of boron, calcium and titanium.

We continue the practice, pioneered by us several years ago, of providing with all Spectrographically Standardised Substances analytical reports stating the quantity of each impurity detected in the particular batch.

Johnson  Matthey

JOHNSON, MATTHEY & CO., LIMITED
3-83 HATTON GARDEN, LONDON, E.C.1
Telephone: Holborn 6989





CONSTANT LOAD WEIGHING BY SUBSTITUTION



As an alternative to the applied-load system of weighing, Stanton now offers the Unimatic single-pan balances—constructed on the constant-load principle of weighing by substitution. The compact, sturdy aluminium case enclosing the instrument is finished in green and black hammertone. Width 13", depth 16", height 20½", weight 32 lb. With the introduction of the Unimatic series, the Stanton range of precision balances now provides an unrivalled choice to meet every requirement.

- | | |
|--------------------------|---|
| UNIMATIC Balances | maintain their accuracy over the whole weighing range being designed on the now well known principle of weighing by substitution. |
| UNIMATIC Balances | are simple to use—even an inexperienced operator can learn how in a few minutes. |
| UNIMATIC Balances | are safe in use—their non-magnetic stainless steel weights are adjusted to within the new N.P.L. tolerances for analytical weights (errors less than 1 in 200,000). They are well protected within the balance case and their accuracy cannot vary. |
| UNIMATIC Balances | are fast—very fast in use. Their controls are conveniently positioned at bench level. Time and Motion study experts have approved this design which results in less fatigue and more weighings in less time—in fact less than 30 seconds for the weighing of an unknown object. |

STANTON UNIMATIC BALANCES

Write for details to STANTON INSTRUMENTS LTD. 119 Oxford Street, London W.1

Telephone: Gerrard 7533 Cables: Stanbal, London, W.1



SOUTHERN-HARWELL PULSE POLAROGRAPH

Following his work on square-wave polarography, Dr. G. C. Barker, of Harwell, has developed a pulse technique which gives greater sensitivity to irreversibly reduced ions, greater independence of base electrolyte impurities, and polarograms of greater certainty and ease of measurement.

The Southern-Harwell Pulse Polarograph provides both normal and derivative polarograms, and is the most sensitive polarograph available.

SOUTHERN ANALYTICAL LIMITED

FRIMLEY ROAD
CAMBERLEY
SURREY
ENGLAND



ONE SCAN SWITCH TWO SPEEDS THREE MINUTES

INFRACORD

Model 137 Spectrophotometer

Versatile, easy to use
Only three controls
Manual setting of drum
No run-back time
Wavelength range: 2.5—15 μ
(Also available for 12.5—25 μ)
Reproducibility: 0.01 μ
Double-beam, aspherical optics
Fully automatic operation
Full range of accessories available
Size: 28" x 18" x 16½"
(71 x 46 x 41 cm)
Made in Great Britain

Perkin-Elmer's famous Infracord Spectrophotometers now being produced at Beaconsfield include yet two further refinements, giving this most-used of all infrared instruments a still greater lead as the world's most practical and economical spectrophotometer: The single scan switch allows choice of *two* speeds, 3 minutes for rapid survey work and 12 minutes for quantitative runs. Also, a wide range of slit programmes are available for greater versatility. The Infracord Spectrophotometer is quick and easy in use, and fully automatic in operation. It takes up little room, and the price is most reasonable. The ideal instrument to provide accurate and reliable information in the widest range of applications in research, analysis, production, control and teaching. We should like to send you full details.

Perkin-Elmer

ANALYTICAL INSTRUMENTS

Beaconsfield, Bucks, England. Telephone: Beaconsfield 2020



A most distinguished range of potentiometric recorders



—for ex-stock delivery

It is a strange fact that however early an installation is planned, when it comes to the recorders, they are often wanted in a hurry. AEI 10" recorders are available ex-stock. Equally important is their specification and performance which makes them ideal for all forms of qualitative and quantitative analysis in general scientific applications and Publication No. 2443-1 gives full information which is important to all users and potential users of potentiometric recorders.

For further information contact your local AEI district office or:—

AEI Instrumentation Division, Process Control Department, P.O. Box One, Harlow, Essex



**Instrumentation Division
Associated Electrical Industries Limited**

ALSO: D.C. AMPLIFIERS

REFERENCE UNITS




COLD JUNCTION UNITS

PULSE GENERATORS

SC/145C



PHILIPS TEMPERATURE ATTACHMENTS FOR DIFFRACTOMETRY



Philips range of X-ray Analytical equipment is increased by the addition of two new goniometer attachments; for both high and low temperature diffraction studies.

LOW TEMPERATURE ATTACHMENT PW.1076

For use with Philips wide-range goniometer in the study of powder specimens at temperatures in the 80°K range.

HIGH/LOW TEMPERATURE ATTACHMENT X86N

For investigation of powder, wire or sheet specimens under conditions of high vacuum (10⁻⁵ mm. or better), gas atmosphere, inert atmosphere or small positive pressures, the new furnace will operate at temperatures from ambient to approx. 2100°C., or with liquid nitrogen coolant down to 80°K.

For full details write today to sole distributors in the U.K.:

RESEARCH & CONTROL INSTRUMENTS LTD.

Instrument House, 207 King's Cross Rd., London, W.C.1 Telephone TERminus: 2877

(RCL0113A)



** Quicker than an analytical balance
more accurate than a simple balance*

TORBAL

PRECISION BALANCES

- * NO KNIFE-EDGES * NO BEARINGS * NO SERVICING
- * EASY DIAL OPERATION * ROBUST * EASILY PORTABLE

We are now demonstrating these remarkable American balances all over the country. They give you simple, rapid weighing with direct reading and being very robust can be moved about like no other balance, without loss of accuracy. The complete range includes Model DLT₅, capacity 500g., accuracy 50 mg., price £61. 0. 0; and DLT₂, capacity 120 g., accuracy 5 mg., price £67. 10. 0.

Write now for full details and a demonstration on your bench.



BAIRD & TATLOCK (LONDON) LIMITED, CHADWELL HEATH, ESSEX, ENGLAND
Branches in London, Manchester, Glasgow.

Agents in U.K. and all over the world

TAS/BT.86

AN INVITATION TO ADVANCE

Anderman & Company Ltd.

CORDIALLY INVITE YOU TO VISIT THEM AT STAND No. 13 AT THE FEIGL ANNIVERSARY SYMPOSIUM, BIRMINGHAM UNIVERSITY ON APRIL 10th TO 13th, AND SEE OUR DISPLAY OF:—

- **Royal Berlin Porcelain** — the world's finest porcelain. Originally made to the exacting demands of kings and still satisfying the more exacting demands of modern science.
- **Merck Chemicals and Indicators** — A range of 6,000 chemicals of extra purity for various purposes. Special display for chromatography.
- **Brand Volumetric Ware** — with all its added advantages including acid proof graduation. Brand ware enables you to bring to your work that extra degree of care and accuracy so necessary in this highly scientific age.
- **The Anderman High Torque Stirrer** — with the best performance at the least expense.
- **The Anderman Cutting Bench** — an entirely new device for cutting glass and ceramic tubing etc.
- **Refractory Tubing and Shapes** — a wide range of sizes and materials.
- **Borosilicate Glassware of Quality.**

IF YOU ARE UNABLE
TO ATTEND WRITE
FOR FURTHER
DETAILS TO

**ANDERMAN
&
COMPANY LTD.**

BATTLEBRIDGE HOUSE
87/95 TOOLEY STREET
LONDON, S.E.1

Tel. HOP 1169
0035/6

micro Hammer Mill

for rapid
analytical
grinding



This compact mill will pulverise samples of a few grammes with practically complete recovery for analytical work. Only 30 to 120 seconds required for disintegration of majority of substances. Crushes stone, coke, bone, chemicals, fertilizers, wood, grain, wool, cotton, roots, hay, straw, etc. Stainless steel grinding mechanism. Variable speed motor for better grinding control. Quickly interchangeable slide-in screens for accurate particle size adjustment. Dust-proof closure. Mill swings wide open with a mere turn of a thumb-nut allowing thorough cleaning. Maximum charge is 120 cc.

Also available:

Micro Mixer Mill for very rapid ultra-fine grinding, mixing, and homogenizing of samples under 1 gramme, for analytical chemistry, spectroscopy, X-ray diffraction, etc. Larger Mixer Mill with

similar applications for samples up to 20 gm Analytical Micro Pulverizer with cooling jacket. Other small scale laboratory mills.

GLENCRESTON LTD.

The Redhouse, 37 Broadway, Stanmore, Middlesex. Telephone: Grimsdyke 4218